



FLAT PANEL DIGITAL RADIOGRAPHY

Let Raymax imaging equipment show you what lies beyond...



Portable Digital X-Ray Panel
Varian 14 x 17

Fixed Digital X-Ray Panel
Toshiba 17 x 17

Raymax is an ISO9001 & ISO13485
Registered Company

Rotate/Transverse Tubestand

- Floor-to-wall or floor-to-ceiling mount with 8 ft. rails
- Transverse tubearm +/- 5", platform tube mount
- Column rotation
- Angulation dial, operator handgrips, and electric locks
- Operates from 24 VDC, 60 Hz, 2 Amp electrical service



Elevator Table

- Four-way float-top table with foot treadle control
- +/- 5 transverse ravel, 32" longitudinal travel
- 84" x 30" table top, variable table height, 22" to 33"
- Foot treadle operation of table top and height movements
- Grid cabinet with electric lock, 22" ravel
- Grid: 8:1 ratio, 103 LPI, Heavy-duty cassette tray
- 400lb. rated patient load
- Operates from 120 VAC, 60 Hz, 6 Amp electrical service

Wallstand Package

- Floor-to-wall mounted column
- Grid cabinet with electric lock
- Grid: 8:1 ratio, 103 LPI, Heavy duty cassette tray
- External 24 VDC, 60 Hz, 2 Amp power supply (120 VAC input)



*Appearance and specifications subject to change without notice.

DISTRIBUTED BY:

Raymax Medical Corp.
Brampton, Ontario, Canada L6T 4P7
Toll Free: (866) RAYMAX1
Tel.: (905) 791 3020
Fax: (905) 791 3375
Email: sales@raymaxmedical.com
www.raymaxmedical.com