

## STORED ENERGY GENERATOR – SE 30

- NO BATTERIES!
- 30KW FROM A 240 VAC, 15 AMP RECEPTACLE

Utilizing the latest in high density electrolytic capacitor storage techniques, the SE-30 has no batteries to slowly lose their charge and require expensive replacement. Rated at more than 10,000 hours of continuous usage, these space age capacitors are at the heart of the stored energy solution provided by the SE-30.



**Stored Energy High Frequency X-ray Generator** offers a “no compromise” solution for customers who have a power problem. The 30 kW SE-30 will provide all the x-ray energy needed for your Veterinary, Chiropractic or Medical Clinics. Plug it in to a 240 VAC, 15 amp outlet and you’re ready to calibrate and take x-rays.

20 Strathearn Avenue Unit 3  
Brampton, Ontario, Canada L6T 4P7

[www.raymaxmedical.com](http://www.raymaxmedical.com)  
[sales@raymaxmedical.com](mailto:sales@raymaxmedical.com)

Toll Free: 1-866-RAYMAX1  
Tel: 1-905-791-3020  
Fax: 1-905-791-3375



### Technical Information:

Output power – 30kW  
Technique factors – 400 mA, 125 kVp, 30 mAs (TV) 100 mAs (TH)  
mA Range – 25 to 400 mA  
kVp Range – 40 to 125 kVp  
Time – 2ms to 2 sec.  
Input power – 240 VAC, Single Phase, 15 Amp Breaker  
Size – 44" x 24" x 12.5"  
Weight – 450 lbs (205 kgs)

### Features:

APR – optional  
Bucky Control – optional  
Tube Protection – Yes  
Note: If the exact mAs is not selectable, choose the next closest.



20 Strathearn Avenue Unit 3  
Brampton, Ontario, Canada L6T 4P7

[www.raymaxmedical.com](http://www.raymaxmedical.com)  
[sales@raymaxmedical.com](mailto:sales@raymaxmedical.com)

Toll Free: 1-866-RAYMAX1  
Tel: 1-905-791-3020  
Fax: 1-905-791-3375