




RAYMAX

MEDICAL IMAGING



Chiropractic Direct Digital X-Ray



**Let Raymax imaging equipment
show you what lies beyond...**

Raymax equipment meets all North American Regulatory Standards



RAYMAX

MEDICAL IMAGING

- Operating Modes: kV/mAs, kV/mA/Time
- 32/40/50/65/80 & 100 kW power levels
- Compatible with both digital and conventional radiographic (RAD) systems
- kV/mAs Stepping



Generator:

Generator Type

CPI 40 kW
High Frequency
Output up to 400 kHz
40 - 150 kV
10 - 500 mA
1 Phase, 230 VAC
768 Techniques
kV/mAs

kV (Radiographic)
mA Range RAD
Line Nominal/Phase²
Anatomical Programs RAD
Technique Selection

Approvals:

CSA, U/L, FDA, Health Canada

DISTRIBUTED BY:



RAYMAX

MEDICAL IMAGING

Brampton, Ontario, Canada L6T 4P7
Toll Free: (866) RAYMAX1
T: (905) 791-3020
F: (905) 791-3375
Email: sales@raymaxmedical.com
www.raymaxmedical.com